

5-Acetylsalicylic acid

Catalog No: #S2735


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	5-Acetylsalicylic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	180.16
Formulation	C ₉ H ₈ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC(=O)C1=CC(C(O)=O)=C(O)C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ulu N, et al. Turk J Med Sci. 2015;45(4):812-9.

Note: This product is for in vitro research use only