

Sodium Demethylcantharidate

Catalog No: #S2736


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sodium Demethylcantharidate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	230.14
Formulation	C ₈ H ₈ Na ₂ O ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C₁CC2C(C(C1O2)C(=O)[O-])C(=O)[O-].[Na+].[Na+]
Pathways: Others
Receptor: Others
Boiling pt: Melting
pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Goldstein P, Leshem M. *Appetite*. 2014 Aug;79:83-90

Note: This product is for in vitro research use only