

Neohesperidin

Catalog No: #S2737


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neohesperidin
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others inhibitor
Calculated MW	610.56
Formulation	C28H34O15
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)O)[C@@H]2[C@H]([C@@H]([C@H](OC2OC3=CC(=C4C(=O)CC(OC4=C3)C5=C(C(=C(C=C5)OC)O)O)CO)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 122 mg/mL (199.81 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lee JH, et al. PhytOthers Res, 2009, 23(12), 1748-1753.

Note: This product is for in vitro research use only