

## Dehydroandrographolide

Catalog No: #S2738


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Dehydroandrographolide
Brief Description	Inhibitors
Purification	99.93%
Target Name	NOS activator
Calculated MW	348.43
Formulation	C <sub>20</sub> H <sub>28</sub> O <sub>5</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: CC1(CO)C(O)CCC2(C)C(CC=C3C=COC3=O)C(=C)CCC12  
 Pathways: Immunology/Inflammation  
 Receptor: iNOS  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
 Appearance: White Needle Crystal  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Hsieh MJ, et al. Oncotarget. 2015 Oct 13;6(31):30831-49.

Note: This product is for in vitro research use only