

Shikimic Acid

Catalog No: #S2739


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Shikimic Acid
Brief Description	Inhibitors
Purification	99.39%
Target Name	Others
Calculated MW	174.15
Formulation	C7H10O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:C1[C@H]([C@@H]([C@@H](C=C1C(=O)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Water: 35 mg/mL (200.97 mM); DMSO: 35 mg/mL (200.97 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Brown SA, et al. Nature, 1955, 175(4459), 688-689.

Note: This product is for in vitro research use only