

## Glycyrrhethinic acid(Enoxolone)

Catalog No: #S2745


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Glycyrrhethinic acid(Enoxolone)
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	470.68
Formulation	C30H46O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:C[C@]12CC[C@](C[C@@H]1C3=CC(=O)[C@@H]4[C@]5(CC[C@@H](C[C@@H]5CC[C@]4([C@@]3(CC2)C)C)(C)O)C)(C)C(=O)O

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Parida PK, et al. Bioorg Med Chem Lett. 2014 Aug 15;24(16):3865-8.

Note: This product is for in vitro research use only