

## Endothall

Catalog No: #S2747


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Endothall
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	186.16
Formulation	C8H10O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:OC(=O)C1C2CCC(O2)C1C(O)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Boknk P, et al. Pharmacology. 2000 Jul;61(1):43-50.

Note: This product is for in vitro research use only