

HIF-1 α Polyclonal Antibody

Catalog No: #41005

Package Size: #41005-1 50ul #41005-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HIF-1 α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF/ICC
Species Reactivity	Hu Ms Rt
Specificity	HIF-1 α Polyclonal Antibody detects endogenous levels of HIF-1 α protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human HIF-1 α .
Target Name	HIF-1 α
Other Names	HIF1A; BHLHE78; MOP1; PASD8; Hypoxia-inducible factor 1-alpha; HIF-1-alpha; HIF1-alpha; ARNT-interacting protein; Basic-helix-loop-helix-PAS protein MOP1; Class E basic helix-loop-helix protein 78; bHLHe78; Member of PAS protein 1; PAS doma
Accession No.	Swiss-Prot: Q16665NCBI Gene ID: 3091
Uniprot	Q16665
GeneID	3091;
SDS-PAGE MW	90kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of HELA cells using HIF-1α Polyclonal Antibody

Note: This product is for in vitro research use only