

Ginkgolide A

Catalog No: #S2749


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ginkgolide A
Brief Description	Inhibitors
Purification	98.00%
Target Name	GABA Receptor antagonist
Calculated MW	408.4
Formulation	C ₂₀ H ₂₄ O ₉
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: CC1C(=O)OC2CC34C5CC(C(C)(C)C)C33C(O)C(=O)OC3OC4(C(=O)O5)C12O
 Pathways: Membrane transporter/Ion channel; Neuroscience
 Receptor: GABAR
 Boiling pt: Melting pt: Solubility: Ethanol: 2 mg/mL (4.89 mM); DMSO: 81 mg/mL (198.33 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhao Q, Gao C, Cui Z. Int Immunopharmacol. 2015 Apr;25(2):242-8.

Note: This product is for in vitro research use only