

## Ginkgolide C

Catalog No: #S2750


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ginkgolide C
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	440.4
Formulation	C <sub>20</sub> H <sub>24</sub> O <sub>11</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other  
 SMILES: CC1C(=O)OC2C(O)C34C5OC(=O)C3(OC3OC(=O)C(O)C43C(C5O)C(C)(C)C)C12O  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

- Liou CJ, et al. Evid Based Complement Alternat Med. 2015;2015:298635.

Note: This product is for in vitro research use only