

Ginkgolide B

Catalog No: #S2751


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ginkgolide B
Brief Description	Inhibitors
Purification	98.00%
Target Name	PAFR antagonist
Calculated MW	424.4
Formulation	C ₂₀ H ₂₄ O ₁₀
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: CC1C(=O)OC2C(O)C34C5CC(C(C)(C)C)C33C(O)C(=O)OC3OC4(C(=O)O5)C12O
 Pathways: GPCR/G
 Protein Receptor: PAFR
 Boiling pt: Melting pt: Solubility: DMSO: 85 mg/mL (200.28 mM)
 Appearance: White Crystal
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Liu X, et al. Thromb Res. 2014 Nov;134(5):1066-73.

Note: This product is for in vitro research use only