

2'-Deoxyadenosine monohydrate

Catalog No: #S2752


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	2'-Deoxyadenosine monohydrate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	269.26
Formulation	C10H15N5O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:C1[C@@H]([C@H](O[C@H]1N2C=NC3=C2N=CN=C3N)CO)O.OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kore AR, et al. Nucleosides Nucleotides Nucleic Acids. 2015;34(1):33-9.

Note: This product is for in vitro research use only