

Loganin

Catalog No: #S2759


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Loganin |
| Brief Description | Inhibitors |
| Purification | 99%+ |
| Target Name | NF-κB inhibitor |
| Calculated MW | 390.38 |
| Formulation | C17H26O10 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research

Area: Immunology SMILES: C[C@H]1[C@H](C[C@H]2[C@@H]1[C@@H](OC=C2(=O)OC)O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)O Pathways: NF-κB Receptor: NF-κB Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Liu K, et al. Life Sci. 2015 Feb 15;123:78-85.

Note: This product is for in vitro research use only