

Naringin dihydrochalcone

Catalog No: #S2761


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Naringin dihydrochalcone
Brief Description	Inhibitors
Purification	99.16%
Target Name	Others
Calculated MW	582.55
Formulation	C27H34O14
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)O[C@@H]2[C@H]([C@@H]([C@H](O[C@H]2OC3=CC(=C(C(=C3)O)C(=O)CCC4=CC=C(C(=C4)O)CO)O)O)O)OPathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Ethanol: 100 mg/mL (171.65 mM); DMSO: 100 mg/mL (171.65 mM) Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Shin W, et al. J Med Chem, 1995, 38(21), 4325-4331.

Note: This product is for in vitro research use only