

Kasugamycin hydrochloride

Catalog No: #S2762


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Kasugamycin hydrochloride
Brief Description	Inhibitors
Purification	90.00%
Target Name	Antifungal agonist
Calculated MW	415.83
Formulation	C ₁₄ H ₂₆ CIN ₃ O ₉
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C[C@@H]1[C@H](C[C@@H]([C@H](O1)OC2[C@@H]([C@H](C[C@@H]([C@@H]2O)O)O)O)N)N=C(C(=O)O)N. CIPat
 hways: Microbiology & virology Receptor: Antifungal Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a
 higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several
 month.

References

- Schluenzen F, et al. Nat. Struct. Mol. Biol. 13, 871C878

Note: This product is for in vitro research use only