

Rosmarinic acid

Catalog No: #S2765


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rosmarinic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	MAO inhibitor; Transferase inhibitor
Calculated MW	360.33
Formulation	C18H16O8
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: C1=CC(=C(C=C1C[C@H](C(=O)O)OC(=O)/C=C/C2=CC(=C(C=C2)O)O)O
 Pathways: Metabolism
 Receptor: MAO-A; MAO-B; COMT
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lee JW, et al. Mol Nutr Food Res. 2015 Mar;59(3):386-400.

Note: This product is for in vitro research use only