

## Neohesperidin Dihydrochalcone (Nhdc)

Catalog No: #S2767


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Neohesperidin Dihydrochalcone (Nhdc)
Brief Description	Inhibitors
Purification	99.80%
Target Name	Others
Calculated MW	612.58
Formulation	C28H36O15
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)O[C@@H]2[C@H]([C@@H]([C@H](O[C@H]2OC3=CC(=C(C(=C3)O)C(=O)CCC4=CC(=C(C=C4)OC)O)O)CO)O)O)OP Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Ethanol: 100 mg/mL (163.24 mM); DMSO: 100 mg/mL (163.24 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Kroeze JH, Chem Senses, 2000, 25(5), 555-559.

Note: This product is for in vitro research use only