

Saikosaponin A

Catalog No: #S2768


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Saikosaponin A
Brief Description	Inhibitors
Purification	98.00%
Target Name	NF-κB inhibitor
Calculated MW	780.99
Formulation	C42H68O13
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C[C@@H]1[C@@H]([C@@H]([C@H]([C@@H](O1)O[C@H]2CC[C@]3([C@H]([C@]2(C)CO)CC[C@@]4([C@@H]3C=C[C@@]56[C@]4(C[C@@H]([C@@]7([C@H]5CC(CC7)(C)C)CO6)O)C)C)O)[C@H]8[C@@H]([C@H]([C@@H]([C@H](O8)CO)O)O)O

Pathways: NF-κB Receptor: NF-κB

Boling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhao H, et al. Exp Ther Med. 2015 Aug;10(2):823-827.

Note: This product is for in vitro research use only