

## Histamine H1 Receptor Polyclonal Antibody

Catalog No: #41007

Package Size: #41007-1 50ul #41007-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Histamine H1 Receptor Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Human,Mouse,Rat
Specificity	Histamine H1 Receptor Polyclonal Antibody detects endogenous levels of Histamine H1 Receptor protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human Histamine H1 Receptor.
Target Name	Histamine H1 Receptor
Other Names	HRH1; Histamine H1 receptor; H1R; HH1R
Accession No.	Swiss-Prot: P35367NCBI Gene ID: 3269
Uniprot	P35367
GeneID	3269;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of COLO205 cells using Histamine H1 Receptor Polyclonal Antibody

---

Note: This product is for in vitro research use only