

## 3-Aminopropionitrile fumarate

Catalog No: #S2769


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	3-Aminopropionitrile fumarate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	114.12
Formulation	C3H6N1/2C4H4O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:NCCC#N.OC(=O)C=CC(O)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Jia LX, et al. J Pathol. 2015 Jul;236(3):373-83.

Note: This product is for in vitro research use only