

Gentiopiricin

Catalog No: #S2770


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Gentiopiricin
Brief Description	Inhibitors
Purification	99.54%
Target Name	Others
Calculated MW	356.33
Formulation	C16H20O9
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C=C[C@H]1[C@@H](OC=C2C1=CCOC2=O)O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Zhang K, et al. Int J Nanomedicine. 2013;8:3227-39.

Note: This product is for in vitro research use only