

Orcinol glucoside

Catalog No: #S2771


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Orcinol glucoside
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	286.28
Formulation	C13H18O7
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CC1=CC(=CC(=C1)OC2[C@@H]([C@H]([C@@H]([C@H](O2)CO)O)O)O
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ge JF, et al. Eur Neuropsychopharmacol. 2014 Jan;24(1):172-80.

Note: This product is for in vitro research use only