

S-allyl-L-cysteine


Catalog No: #S2772

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	S-allyl-L-cysteine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	161.22
Formulation	C6H11NO2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:C=CCSC[C@@H](C(=O)O)NPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Water:

>5mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Xu X, et al. J Food Sci Technol. 2015 Mar;52(3):1724-9.

Note: This product is for in vitro research use only