

Lycorine chloride

Catalog No: #S2774


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lycorine chloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	323.17
Formulation	C16H18ClNO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C1CN2CC3=CC4=C(C=C3[C@H]5[C@H]2C1=C[C@@H]([C@H]5O)O)OCO4.Cl Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Hu M, et al. Oncotarget. 2015 Jun 20;6(17):15348-61.

Note: This product is for in vitro research use only