

(+)-Catechin Hydrate


Catalog No: #S2778

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	(+)-Catechin Hydrate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Antioxidant agonist
Calculated MW	308.28
Formulation	C15H16O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C1[C@@H]([C@H](OC2=CC(=CC(=C21)O)O)C3=CC(=C(C=C3)O)O)O Pathways: oxidation-reduction Receptor: Antioxidant
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Mehra P, et al. Indian J Exp Biol. 2013 Oct;51(10):823-7.

Note: This product is for in vitro research use only