

Olaquinox

Catalog No: #S2779


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Olaquinox
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	263.25
Formulation	C ₁₂ H ₁₃ N ₃ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CC1=C(C(=O)NCCO)[N+](=O)C2=CC=CC=C2N1[O-]
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting
 pt: Solubility: DMSO: Soluble
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhao WX, et al. Cell Biol Toxicol. 2013 Aug;29(4):229-38.

Note: This product is for in vitro research use only