

Catalpol

Catalog No: #S2780


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Catalpol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	362.33
Formulation	C15H22O10
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C1=CO[C@H]([C@H]2[C@@H]1[C@@H]([C@H]3[C@@]2(O3)CO)O)[C@H]4[C@@H]([C@H]([C@H]([C@H](O4)CO)O)O)O

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang JH, et al. Biomed Pharmacol. 2015 Feb;69:291-6.

Note: This product is for in vitro research use only