

## Sanguinarine

Catalog No: #S2781


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sanguinarine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Antioxidant agonist
Calculated MW	332.33
Formulation	C <sub>20</sub> H <sub>14</sub> NO <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: C[N+]1=C2C3=CC4=C(OCO4)C=C3C=CC2=C2C=CC3=C(OCO3)C2=C1 Pathways: oxidation-reduction Receptor: Antioxidant  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Mehra P, et al. Indian J Exp Biol. 2013 Oct;51(10):823-7.

Note: This product is for in vitro research use only