

Methyl?Dihydrojasmonate;Hedione

Catalog No: #S2783


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Methyl?Dihydrojasmonate;Hedione
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	226.32
Formulation	C13H22O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:CCCCC1C(CC(=O)OC)CCC1=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. da Silva GB, et al. Int J Nanomedicine. 2014 Feb 11;9:867-76.

Note: This product is for in vitro research use only