

Orsellinic acid ethyl ester

Catalog No: #S2784


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Orsellinic acid ethyl ester
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	196.2
Formulation	C10H12O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCOC(=O)C1=C(O)C=C(O)C=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Qi H, Chen J, Shi Y. Zhongguo Zhong Yao Za Zhi. 2011 Aug;36(16):2233-5. Chinese.

Note: This product is for in vitro research use only