

## Oxysophocarpine

Catalog No: #S2786

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Oxysophocarpine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Prostaglandin Receptor agonist
Calculated MW	262.35
Formulation	C15H22N2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: C1C[C@H]2CN3[C@H](CC=CC3=O)[C@@H]4[C@H]2[N@@+](C1)(CCC4)[O-] Pathways: GPCR/G  
 Protein Receptor: PGE2 Binding pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Yang Y, et al. *Planta Med.* 2015 Jul;81(10):791-7.

Note: This product is for in vitro research use only