

Picroside I

Catalog No: #S2787


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Picroside I
Brief Description	Inhibitors
Purification	99%+
Target Name	Others
Calculated MW	492.47
Formulation	C24H28O11
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: OCC1OC(OC2OC=CC3C(OC(=O)C=CC4=CC=CC=C4)C4OC4(CO)C23)C(O)C(O)C1O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Zhu J, et al. J Chromatogr B Analyt Technol Biomed Life Sci. 2015 Jul 1;993-994:47-59.

Note: This product is for in vitro research use only