

## Baccatin III

Catalog No: #S2788


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Baccatin III
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	586.64
Formulation	C31H38O11
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:CC1=C2[C@H](C(=O)[C@@]3([C@H](C[C@@H]4[C@]([C@H]3[C@@H]([C@@](C2(C)C)(C[C@@H]1O)O)OC(=O)C5=CC=CC=C5)(CO4)OC(=O)C)O)C)OC(=O)C

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White

PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Han F, et al. J Sci Food Agric, 2014, 94: 2376C2383.

Note: This product is for in vitro research use only