

## Scutellarin

Catalog No: #S2789


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Scutellarin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	462.37
Formulation	C21H18O12
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:C1=CC(=CC=C1C2=CC(=O)C3=C(C(=C(C=C3O2)O[C@@H]4[C@@H]([C@H]([C@@H]([C@H](O4)C(=O)O)O)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Faint Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Wang WW, et al. J Mol Neurosci. 2015 Oct 29.

Note: This product is for in vitro research use only