

Pectolarin

Catalog No: #S2791


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pectolarin
Brief Description	Inhibitors
Purification	99%+
Target Name	Others
Calculated MW	622.59
Formulation	C29H34O15
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)OC[C@H]2[C@H]([C@@H]([C@H]([C@@H](O2)OC3=C(C(=C4C(=C3)OC(=C4=O)C5=CC=C(C=C5)OC)O)OC)O)O)O)O)O)O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Yoo YM, et al. Biol Pharm Bull. 2008 Apr;31(4):760-4.

Note: This product is for in vitro research use only