

## Tetrahydropalmatine

Catalog No: #S2793


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Tetrahydropalmatine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	355.44
Formulation	C <sub>21</sub> H <sub>25</sub> NO <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: COC1=C(OC)C=C2C3CC4=C(CN3CCC2=C1)C(OC)=C(OC)C=C4  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt:  
 Solubility: DMSO: 10 mM  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. XING Shu-Hua, et al. Acta Pharmacologica Sinica, 1994-01

Note: This product is for in vitro research use only