

Esculetin


Catalog No: #S2798

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Esculetin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	178.14
Formulation	C ₉ H ₆ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC1=C(O)C=C2C=CC(=O)OC2=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Faint Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Cho JH, et al. Int J Oncol. 2015 Jan;46(1):265-71.

Note: This product is for in vitro research use only