

(20S)-Protopanaxadiol

Catalog No: #S2799


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	(20S)-Protopanaxadiol
Brief Description	Inhibitors
Purification	95.00%
Target Name	Apoptosis inducer inhibitor
Calculated MW	460.39
Formulation	C30H52O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer
 SMILES: CC(=CCC[C@@](C)([C@H]1CC[C@@]2([C@@H]1[C@@H](C[C@H]3[C@]2(CC[C@@H]4[C@@]3(CC[C@@H](C4(C)C)O)C)O)C)O)C
 Pathways: Apoptosis Receptor; Apoptosis inducer
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Needle

Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Liu R, et al. *Fitoterapia*. 2010 Dec;81(8):1079-87.

Note: This product is for in vitro research use only