

## L-Theanine


Catalog No: #S2800

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	L-Theanine
Brief Description	Inhibitors
Purification	99.81%
Target Name	Others
Calculated MW	174.2
Formulation	C7H14N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CCNC(=O)CC[C@@H](C(=O)O)N  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Prez-Vargas JE, et al. Hum Exp Toxicol. 2015 Apr 7. pii: 0960327115578864.

Note: This product is for in vitro research use only