

## Panaxtriol

Catalog No: #S2804


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Panaxtriol
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	476.39
Formulation	C30H52O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:C[C@@]1(CCCC(O1)(C)C)[C@H]2CC[C@@]3([C@@H]2[C@@H](C[C@H]4[C@]3(C[C@H]([C@@H]5[C@@]4(CC[C@@H](C5(C)C)O)C)O)C)O)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Cheng Y, Shen LH, Zhang JT. Acta Pharmacol Sin. 2005 Feb;26(2):143-9.

Note: This product is for in vitro research use only