

Chlorogenic Acid

Catalog No: #S2805


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chlorogenic Acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	354.32
Formulation	C16H18O9
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C1[C@H]([C@H]([C@@H](C[C@@]1(C(=O)O)O)OC(=O)/C=C/C2=CC(=C(C=C2)O)O)O)Pathways:OthersReceptor:OthersBoli
ng pt:Melting pt:Solubility:Ethanol: 71 mg/mL (200.38 mM); Water: 18 mg/mL (50.8 mM); DMSO: 71 mg/mL (200.38 mM)Appearance:White
PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be
stored below -20 °C for several month.

References

1. Illing EA, et al. Otolaryngol Head Neck Surg. 2015 Aug;153(2):291-7.

Note: This product is for in vitro research use only