

Lappaconitine

Catalog No: #S2806


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lappaconitine
Brief Description	Inhibitors
Purification	92.57%
Target Name	Others
Calculated MW	584.64
Formulation	C32H44N2O8
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CCN1C[C@@]2(CC[C@@H]([C@@]34[C@@H]2C[C@H](C31)[C@]5(C[C@@H]([C@H]6C[C@@H]4[C@@]5([C@H]6OC)O)C)O)OC(=O)C7=CC=CC=C7NC(=O)C Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Bryzgalov AO, et al. Cardiovasc Hematol Agents Med Chem. 2013 Sep;11(3):211-7.

Note: This product is for in vitro research use only