

## Caffeic Acid

Catalog No: #S2807


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Caffeic Acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	180.16
Formulation	C <sub>9</sub> H <sub>8</sub> O <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: C1=CC(=C(C=C1/C=C/C(=O)O)O)O  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: DMSO: 36 mg/mL (199.82 mM)  
Appearance: Khaki Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Pradhananga S, Shim WS. Eur J Pharmacol. 2015 Sep 5;762:313-21.

Note: This product is for in vitro research use only