

Adenosine disodium triphosphate

Catalog No: #S2812


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adenosine disodium triphosphate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	551.15
Formulation	C ₁₀ H ₁₄ N ₅ Na ₂ O ₁₃ P ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: C1=NC2=C(C(=N1)N)N=CN2[C@H]3[C@@H]([C@@H]([C@H](O3)COP(=O)(O)OP(=O)([O-])OP(=O)(O)[O-])O)O.[Na+].[Na+]Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Shan G, et al. Se Pu. 2014 Nov;32(11):1275-9. Chinese.

Note: This product is for in vitro research use only