

## Puerarin (Kakonein)

Catalog No: #S2815


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Puerarin (Kakonein)
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others antagonist
Calculated MW	416.38
Formulation	C <sub>21</sub> H <sub>20</sub> O <sub>9</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: C1=CC(=CC=C1C2=C(C(=O)C=CC(=C3[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)O)OPathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 86mg/ml(206.54 mM); Water:

## References

1. Overstreet DH, et al. Pharmacol Biochem Behav, 2003, 75(3), 619-625.

Note: This product is for in vitro research use only