

Picroside II

Catalog No: #S2817


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Picroside II
Brief Description	Inhibitors
Purification	99.68%
Target Name	Others
Calculated MW	512.46
Formulation	C23H28O13
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:COC1=C(C=CC(=C1)C(=O)OC2C3C=COC(C3C4(C2O4)CO)O[C@H]5[C@@H]([C@H]([C@@H]([C@H](O5)CO)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang L, et al. Exp Ther Med. 2015 Apr;9(4):1253-1258.

Note: This product is for in vitro research use only