

## Methyl Vanillate


Catalog No: #S2818

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Methyl Vanillate
Brief Description	Inhibitors
Purification	99.50%
Target Name	Wnt/beta-catenin activator
Calculated MW	182.17
Formulation	C9H10O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer  
SMILES: COC(=O)C1=CC(OC)=C(O)C=C1  
Pathways: Cytoskeletal Signaling; Stem Cells  
Receptor: Wnt  
Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Cha PH, et al. PLoS One. 2014 Jan 22;9(1):e85546.

Note: This product is for in vitro research use only