

Methyl 4-hydroxycinnamate

Catalog No: #S2819


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Methyl 4-hydroxycinnamate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	178.18
Formulation	C10H10O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC(=O)/C=C/C1=CC=C(C=C1)O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Tan EM, et al. J Phys Chem B. 2013 May 2;117(17):4798-805.

Note: This product is for in vitro research use only