

## Crocin

Catalog No: #S2823


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Crocin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	976.96
Formulation	C44H64O24
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:C/C(=C\C=C\C=C\C(=C\C=C\C(=O)O[C@@H]1O[C@@H]([C@H]([C@@H]([C@H]1O)O)O)CO[C@@H]2O[C@@H]([C@H]([C@@H]([C@H]2O)O)O)CO)/C)/C=C/C(=O)O[C@@H]3O[C@@H]([C@H]([C@@H]([C@H]3O)O)O)CO[C@@H]4O[C@@H]([C@H]([C@@H]([C@H]4O)O)O)CO)\CPATHWAYS:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Red PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Sapanidou V, et al. Theriogenology. 2015 Nov;84(8):1273-82.

Note: This product is for in vitro research use only