

Spiculisporic Acid

Catalog No: #S2828


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Spiculisporic Acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	328.41
Formulation	C17H28O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CCCCCCCCC(C(O)=O)C1(CCC(=O)O1)C(O)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang R, et al. Molecules (Basel, Switzerland),2012, 17 (11): 13175C82.

Note: This product is for in vitro research use only